

# #631-0115 Advanced Microscope

## INFORMATION & INSTRUCTIONS

Please read the booklet before using the microscope.

### MODEL SME-F6 Microscope, Monocular, Light. BIOLOGICAL MICROSCOPE

This information manual is intended to familiarize the first time user with their instrument. These microscope are ruggedly built for classroom and laboratory use. The instruments feature fully-coated achromatic optics and international standard DIN objectives and eyepieces. Optically, they are excellent systems with good resolution, centering, and paragon features.

Mechanical movements are smooth and positive.

Mechanical stage, and rack stop serve to minimize damage.


#### SPECIFICATIONS:

- Eyepiece:

| Type | Magnification | Vision field's diameter |
|------|---------------|-------------------------|
| WF   | 10X           | 18mm                    |
| WF   | 16X           | 11mm                    |

- Viewing Angle: Head inclined at 30°, 360°rotatable.
- Illumination: Built-in 220V/110V20W Incandescent lamp (Max power: 20W), and there are an ON-OFF switch & a fuse for prevention of short circuit. With focusing lens of Abbe condenser (N.A 1.25) and iris diaphragm.
- Objective:

| Type                 | Magnification | N.A   | Working Distance |
|----------------------|---------------|-------|------------------|
| Achromatic Objective | 4X            | 0.10  | 37.5mm           |
|                      | 10X           | 0.25  | 6.54mm           |
|                      | 40X(S)        | .0.65 | 0.63mm           |
|                      | 100(S) Oil    | 1.25  | 0.18mm           |



#### SET-UP:

##### 1. GETTING THE MICROSCOPE READY FOR USE

First, familiarize yourself both visually and physically with the mechanical parts of your microscope. Gently operate each part by hand (no tools required) to see how it behaves and what result it produces.

##### 2. MOUNTING THE OBJECTIVE AND EYEPIECE

To put the eyepiece in place simply slide it into the eyepiece tube. The objective lenses screw into the rotating nosepiece. If the four objectives are already mounted in the nosepiece, check to see if they are tightly in place.

To do this, grip the knurled ring on the objectives with your fingers and turn it to the right. Do not tighten them

with tools, they need only to be “finger tight”. With all four objectives mounted in the nosepiece the entire nosepiece can be rotated by hand. As you rotate the nosepiece you will notice that each objective lens has a “click stop” which stops and holds the objective lens being used in its proper position relative to the optical axis of the microscope.

### 3. ILLUMINATION

In order to provide illumination for your new microscope, simply plug the power cord into an electrical socket and switch your microscope “ON”, and you are ready to begin working.

### OPERATION:

Once you have set up your microscope for use, follow the steps below to prepare the specimen to be studied for viewing:

- Place the slide to be viewed onto the stage so that the stage clips hold it firmly in place. Position the slide so that the portion of the specimen that you wish to observe is below the 4X objective.
- Slowly rotate the coarse focus knob until the stage is approximately 1/8” from the objective, being sure not to let the slide come in contact with the objective, for damage to either or both may occur. When using the coarse focus knob, start with the objective as close to the stage as possible and focus from the bottom up. This will help to prevent any contact between the objective and the stage.
- Adjust the iris diaphragm attached to the condenser to control the angle of light illumination the specimen.
- To change the magnification, first draw the objective to the position furthest from the stage. Then turn the objective turret until the desired objective is in place. Return the objective to the point closest to the stage and focus once again backing away from the stage.
- To determine the total magnification, multiply the magnification of the objective with that of the eyepiece. (e.g., 10X objective X 10X eyepiece = 100X magnification).

### MAINTENANCE:

1. Keep the instrument in dry and cool place, and keep it away from acid or alkali.
2. Put the dust cover over it when it is not used.
3. Keep lenses clean, dust on it is moved with soft brush.
4. When you change the lamp or fuse, you must disconnect the power.
5. This instrument is used in 0-40°C circumstance.
6. Class of dustproof and waterproof: 2P20.
7. Fuse spec: 250V, F=0.5A
7. On the instrument showing drawings show: Disconnect Power Before Opening.

---

**Warranty and Parts:** We replace all defective or missing parts free of charge. Additional replacement parts may be ordered toll-free. We accept MasterCard, Visa, checks and School P.O.s. All products warranted to be free from defect for 90 days. Does not apply to accident, misuse or normal wear and tear. Intended for children 13 years of age and up. This item is not a toy. It may contain small parts that can be choking hazards. Adult supervision is required.